



# Impregnated Dressings for Advanced Wound Care



[medsupplyinc.com](http://medsupplyinc.com)

# Honey Gauze Dressing



REF	UNIT OF MEASURE	CONTENTS	SOLD BY	SIZE
HOD4X4	Box (BX)	10 Each	1 Box of 10	10 × 10cm
HOD6X6	Box (BX)	10 Each	1 Box of 10	15 × 15cm

MedSupply Gauze Honey Dressing is an easy to apply dressing impregnated with 100% sterile, medical-grade manuka honey.

MedSupply Gauze Honey Dressing is specially formulated to treat wounds and help them heal faster, while preventing infections.

The Gauze Honey Dressing secures the wound by absorbing or donating moisture to the wound while providing a cooling, soothing relief.

## Features:

Eliminate odors

Osmotic effect – harmful tissues is draw away from the wound bed

Maintains the ideal moist wound healing environment

## Uses:

- + Venous ulcers
- + Diabetic foot ulcers
- + Pressure ulcers
- + Minor burns

## Contraindications:

Honey Gauze Dressings are not indicated for use in third degree burns or to control heavy bleeding. It is not indicated for patients that have a known sensitivity to honey or any other component parts specific to each dressing.

# Hydrogel Gauze Dressings

Hydrogel Gauze dressings are designed to provide a moist wound-healing environment saturated with an amorphous hydrogel. The dressing size of 10×10cm are impregnated with 1.0 ounce of hydrogel and the 15×15cm are impregnated with 2.0 ounce of hydrogel.

Amorphous hydrogels are formulations of water, polymers and other ingredients with no shape, designed to donate moisture to a dry wound and to maintain a moist healing environment. The high moisture content serves to rehydrate wound tissue.



REF	UNIT OF MEASURE	CONTENTS	SOLD BY	SIZE
HIG4X4	Box (BX)	10 Each	1 Box of 10	10 × 10cm
HIG6X6	Box (BX)	10 Each	1 Box of 10	15 × 15cm

## Features:

Provides optimal environment for moist wound healing

Conforms to wound bed

## Uses:

- + Venous leg ulcers
- + Arterial ulcer
- + Pressure ulcer
- + 1st and 2nd Degree Burns
- + Skin graft donor site
- + Abrasions

## Contraindications:

Not for use on third-degree burns or on excessively draining wounds and infected wounds. Patients with a known hypersensitivity to the product itself or to its components.

# Paraffin Gauze Dressing



REF	UNIT OF MEASURE	CONTENTS	SOLD BY	SIZE
PGD4X4	Box (BX)	10 Each	1 Box of 10	10 × 10cm
PGD6X6	Box (BX)	10 Each	1 Box of 10	15 × 15cm

Oil Emulsion Wound Dressing consists of a non-adherent gauze mesh impregnated with white petrolatum in an oil emulsion blend.

This primary dressing permits the flow of exudate without sticking to the wound bed. Conforms to all body contours.

Non-adherent, use beneath a secondary dressing to manage exudate and securement.

## Features:

Non-adhesive

Sterile

Not made with natural rubber latex.

Atraumatic removal

Help maintain a moist wound environment

Minimal maceration and occlusion due to open mesh gauze

## Uses:

- + Abrasions
- + 1st and 2nd degree burns
- + Surgical Laceration

## Contraindications:

Not for use in patients known to have sensitivity to any of the dressing components. Do not use in cavity wounds.

# Petrolatum Gauze Dressing

Petrolatum Gauze is a fine mesh, absorbent gauze impregnated with white petrolatum. It remains moist, nontoxic and nonirritating. It is non-adherent and sterile.

The non-adherent cotton gauze is very comfortable for the patient to wear. The gauze easily clings and conforms to all the contours of the body for an ideal and comfortable fit.



REF	UNIT OF MEASURE	CONTENTS	SOLD BY	SIZE
PTL4X4	Box (BX)	10 Each	1 Box of 10	10 × 10cm
PTL6X6	Box (BX)	10 Each	1 Box of 10	15 × 15cm

## Features:

Non-adhesive  
Abrasions  
Suture Wounds  
Surgical incisions

## Uses:

- + Minor burns
- + Abrasions
- + Suture wounds
- + Surgical incisions

## Contraindications:

Not for use on third-degree burns or on excessively draining wounds and infected wounds. Patients with a known hypersensitivity to the product itself or to its components.

# Xeroform Gauze Dressing



REF	UNIT OF MEASURE	CONTENTS	SOLD BY	SIZE
XRFG1X8	Box (BX)	10 Each	1 Box of 10	2.5 × 20cm
XRFG2X4	Box (BX)	10 Each	1 Box of 10	05 × 05cm
XRFG4X4	Box (BX)	10 Each	1 Box of 10	10 × 10cm
XRFG6X6	Box (BX)	10 Each	1 Box of 10	15 × 15cm

Sterile fine mesh gauze occlusive dressing impregnated with petrolatum and 3% Xeroform (Bismuth Tribromophenate). It is used to treat everything from minimal drainage to moderate and heavy drainage wounds; such as, surgical incisions, lacerations and burns to skin grafts, donor sites, dermal ulcers, and more. They are non-adherent and help maintain a moist wound environment. Clings and conforms well to body contours.

## Features:

Minimally draining wounds

Surgical incisions

Lacerations

Donor sites

First and second degree burns

## Benefits:

- + Non-adherent
- + Help maintain a moist wound environment
- + Occlusive dressing
- + The fine mesh gauze is impregnated with a blend formula of petrolatum and bismuth tribromophenat that provides bacteriostatic protection.

## Contraindications:

Not for use on third-degree burns or on excessively draining wounds and infected wounds. Patients with a known hypersensitivity to the product itself or to its components.

# Amorphous Hydrogel Dressing

MedSupply Amorphous Hydrogel Dressing is useful for rehydrating sloughy or necrotic tissue and enhancing autolytic debridement.

On minimally draining wounds it provides moist wound environment, conducive to wound healing.



REF	UNIT OF MEASURE	CONTENTS	SOLD BY	SIZE
HWD 15G	Box (BX)	1 Each	1 Box of 1	15G

## Features:

Eliminate odors

Osmotic effect – harmful tissues is draw away from the wound bed

Maintains the ideal moist wound healing environment

## Uses:

- + Venous ulcers
- + Diabetic foot ulcers
- + Pressure ulcers
- + Minor burns

## Contraindications:

There are two types of wounds for which a hydrogel dressing is not advised full-thickness burns and moderate- to high-exuding wounds.

**We provide you  
with the product,  
services, and support  
of the highest quality  
while maintaining  
top-level integrity,  
customer satisfaction,  
and ethical standards.**



For more information on our products visit our website or  
contact us today.

MedSupplyInc.com | 954-342-9720  
1241 Stirling Rd. Unit #116  
Dania Beach, FL 33004